

Envisage – Configuring New Media

New Media

From time to time, Egton will distribute new media to your Envisage server. Media might include pictures, videos and presentations.

New media is made available in your system using a patching system from EMIS, which will ensure that the information *only* gets to your envisage server.

Where to find it

Once delivered you only need to access your Envisage Media Scheduler and check the media 'Groups' screen. The new media should show in the 'ungrouped' folder.

You can preview the media by clicking on the item in the Envisage Media Scheduler, then clicking the 'Preview Media' Button.

How to display the media to your patients

Available on every computer with the Envisage Administrator Module installed is the Envisage Media Scheduler Manual. This gives details on how to add the media to a sequence, and how to have that sequence display on the screen.

Associated documentation

None

No part of this document may be sold, hired, reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying, recording and information storage and retrieval systems for any other purpose than the purchaser's use without the express written permission of Egton.

Every effort is made to ensure that your Egton documentation is up to date, but our commitment to constantly improve our software and systems means that there may have been changes since this document was produced.

Egton is a trading name of Egton Medical Information Systems Limited

• tel 0845 125 5525 • fax 0113 203 1647 • mail@egton.net • www.egton.net

Egton • Stafford House • Unit 2 Leathley Road Industrial Estate • Leathley Road • Hunslet • Leeds • LS10 1BG